

How to test using Postman

This guide will show how to interact with the API using the API tool Postman.

Create a new Collection

1. Select "New" > "Collection"


The API will use collection variables to set base URLs and Tokens.

2. Skip Authorization when creating a new collection

The Ship API uses OAuth 2.0. We also delegate permissions to clients based on intended API calls (audience in OAuth).

In the diagram below, in Postman, you cannot set the **audience** parameter.

There is a work around later in this walk through .



The screenshot shows the 'GET NEW ACCESS TOKEN' dialog in Postman. The fields are as follows:

Field	Value
Token Name	Test Client Token
Grant Type	Client Credentials
Access Token URL	https://coll-8-dev.eu.auth0.com/oauth/token
Client ID	[Redacted]
Client Secret	[Redacted]
Scope	e.g. read:org
Client Authentication	Send as Basic Auth header

At the bottom of the dialog is an orange button labeled 'Request Token'.

3. Set Base URL

Enter the base URL for your testing. In this case the URL is set to **https://ship-api-uat.drop2shop.ie/api**.

EDIT COLLECTION

Name: Coll8 Ship API UAT

Description Authorization Pre-request Scripts Tests **Variables**

These variables are specific to this collection and its requests. [Learn more about collection variables.](#)

	VARIABLE	INITIAL VALUE ⓘ	CURRENT VALUE ⓘ	...	Persist All	Reset All
<input checked="" type="checkbox"/>	base_url	https://ship-api-uat.drop2	https://ship-api-uat.drop2shop.ie/api			
<input checked="" type="checkbox"/>	token	eyJhbGciOiJSUzI1NiIsInR5cC	eyJhbGciOiJSUzI1NiIsInR5cC			
	Add a new variable					

6. Test Ping

There is an endpoint for testing you can communicate with the API. It's a Ping endpoint which responds to **POST** requests.

- Set `{{base_url}}` parameter in the url along with `/ping`
- Add Authorization Header set to `Bearer {{token}}`
- Click "Send"

POST Get New Token POST Ping GET Untitled Request No Environment

Ping Examples (0)

POST `{{base_url}}/ping` Send Save

Params Authorization **Headers (1)** Body Pre-request Script Tests Cookies Code

KEY	VALUE	DESCRIPTION	...	Bulk Edit	Presets
<input checked="" type="checkbox"/> Authorization	Bearer <code>{{token}}</code>				
Key	Value	Description			

Body Cookies (1) Headers (8) Test Results Status: 200 OK Time: 209 ms Size: 427 B Save Download

Pretty Raw Preview JSON

```
1 {
2   "message": "OK"
3 }
```